



Protein Standards for Human Cytokine ELISA Strip

Catalog Number EA-1082

(For Research Use Only)


<u>Tube</u>	<u>Protein</u>	<u>Stock Conc.</u>	<u>Volume</u>
#1	Human VEGF	400ng/ml	10ul
#2	Human EGF	400ng/ml	10ul
#3	Human IL-6	400ng/ml	10ul
#4	Human Resistin	400ng/ml	10ul
#5	Human PAI-1	400ng/ml	10ul
#6	Human IL-12	400ng/ml	10ul
#7	Human IL-13	400ng/ml	10ul
#8	Human Eotaxin	400ng/ml	10ul

Preparation of protein standard dilutions

1. Add 200ul of Diluent buffer the wells of the first strip, and add 100ul of Diluent buffer to the wells of the rest strips according to the following table.
2. See Table 1 for appropriate amount of protein standards to add to the first strip according to the following table.
3. Use multi-channel pipette to mix the dilutions and transfer 100ul to the next dilution wells.
4. Repeat the transfer until 5th strip. Trash 100ul from 5th strip after mixing. Do not transfer to the 6th Strip

Note: Substrate incubation time may vary due to different antibodies reactivity. Stronger signals (Strong blue color) could be stopped early after 5 minutes. Weaker signals should be incubated for 10-30 minutes. Always stop protein standards along with samples from the same row at the same time.

Table 1: Dilution of protein standards



<u>Protein</u>	<u>Added to 1st Strip</u>	<u>1st Strip</u>	<u>2nd Strip</u>	<u>3rd Strip</u>	<u>4th Strip</u>	<u>5th Strip</u>	<u>6th Strip</u>
Human VEGF	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Human EGF	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Human IL-6	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Human Resistin	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Human PAI-1	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Human IL-12	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Human IL-13	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Human Eotaxin	4ul	200ul	100ul	100ul	100ul	100ul	100ul
			1:2	1:4	1:8	1:16	Blank

Table 2: Concentrations of Protein Standard Dilutions

<u>Protein</u>		<u>1:2</u>	<u>1:4</u>	<u>1:8</u>	<u>1:16</u>	<u>Blank</u>
Human VEGF	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Human EGF	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Human IL-6	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Human Resistin	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Human PAI-1	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Human IL-12	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Human IL-13	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Human Eotaxin	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank

<u>Protein</u>		<u>1:2</u>	<u>1:4</u>	<u>1:8</u>	<u>1:16</u>	<u>Blank</u>
Rat TNF α	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat VEGF	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat IL-6	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat FGFb	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat IFN γ	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat Leptin	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat MCP-1	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat Rantes	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank

<u>Protein</u>		<u>1:2</u>	<u>1:4</u>	<u>1:8</u>	<u>1:16</u>	<u>Blank</u>
Rat TNF α	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat VEGF	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat IL-6	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat FGFb	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat IFN γ	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat Leptin	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat MCP-1	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank
Rat Rantes	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	0.25ng/ml	Blank

