



Protein Standards for Rat Oxidative Stress ELISA Strip

Catalog Number EA-1502

(For Research Use Only)

<u>Tube</u>	<u>Protein</u>	<u>Stock Conc.</u>	<u>Volume</u>
#1	Rat TNF α	400ng/ml	10ul
#2	Rat TGF β	400ng/ml	10ul
#3	Rat MCP-1	400ng/ml	10ul
#4	Rat IL-1 α	400ng/ml	10ul
#5	Rat IL-1 β	400ng/ml	10ul
#6	Rat IL-6	400ng/ml	10ul
#7	Rat IL-15	400ng/ml	10ul
#8	Rat VEGF	400ng/ml	10ul

Preparation of prote in standard dilutions

1. Add 200ul of Diluent buffer the wells of the first strip, and add 100ul of Diluent buffer to the wells of the rest strips according to the following table.
2. Add appropriate amount of protein standards to the first strip according to the table 1
3. Use multi-channel pipette to mix the dilutions and transfer 100ul to the next dilution wells.
4. Repeat the transfer until 5th strip. Trash 100ul from 5th strip after mixing. Do not transfer to the 6th Strip

Note: Substrate incubation time may vary due to different antibodies reactivity. Stronger signals (Strong blue color) could be stopped early after 5 minutes. Weaker signals should be incubated for 10-30 minutes. Always stop protein standards along with samples from the same row at the same time.

Table 1: Dilution of prote in standards

<u>Protein</u>	<u>Added to 1st Strip</u>	<u>1st Strip</u>	<u>2nd Strip</u>	<u>3rd Strip</u>	<u>4th Strip</u>	<u>5th Strip</u>	<u>6th Strip</u>
Rat TNF α	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Rat TGF β	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Rat MCP-1	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Rat IL-1 α	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Rat IL-1 β	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Rat IL-6	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Rat IL-15	4ul	200ul	100ul	100ul	100ul	100ul	100ul
Rat VEGF	4ul	200ul	100ul	100ul	100ul	100ul	100ul
			1:2	1:4	1:8	1:16	Blank

Table 2: Concentrations of Protein Standard Dilutions

<u>Protein</u>	<u>Added to 1st Strip</u>	<u>1:2</u>	<u>1:4</u>	<u>1:8</u>	<u>1:16</u>	<u>Blank</u>
Rat TNF α	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Rat TGF β	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Rat MCP-1	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Rat IL-1 α	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Rat IL-1 β	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Rat IL-6	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Rat IL-15	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank
Rat VEGF	8ng/ml	4ng/ml	2ng/ml	1ng/ml	0.5ng/ml	Blank