



Signosis


Innovative Plate Assay Solutions

Product Information

pCEBP-Luc Reporter Vector (LR-2022)

Product

pCEBP-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of CEBP in cells.



NheI **CEBP cis element**

TTTCTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGCATTGCGCAATATTGCGCAATATTGCGCAAT

BglII HindIII

AGATCTGGGTATATAATGGAAGCTTGGCATTCCGGTACTGTTGGTAAAGCCACC

Luciferase gene

ATGGAAGACGCCAAAAACATAAAGAAAGGCCCG

Shelf Life

The shelf life of this product is 12 months.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com