



Signosis


Innovative Plate Assay Solutions

Product Information

pCOUP-TF-Luc Reporter Vector (LR-2029)

Product

pCOUP-TF-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of COUP-TF in cells.



```

                                NheI           COUP-TF cis element
TTTTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGC CAGGTCACAGGTCAGACAGGTCACAGGTCAG
  BglII                               HindIII
AAGATCTGGGTATATAATGGAAGCTTGGCATTCCGGTACTGTTGGTAAAGCCACCATGGAAGACGCCAAAA
                                Luciferase
ACATAAAGAAAGGCCCGCGCCATTCTATCCGCTGGAAGATGGAACCGCTGGAGAGCAACTGCATA

```

Shelf Life

The shelf life of this product is 6 months. Use it before expiration.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com