



Signosis


Innovative Plate Assay Solutions

Product Information

pE2F-Luc Reporter Vector (LR-2031)

Product

pE2F-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of E2F in cells.



```

                                NheI                               E2F cis element
TTTTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGCTTTCGCGCTTTTCGCGCTTTTCGCGCTTTTCGCG
  BglII                               HindIII
CAGATCTGGGTATATAATGGAAGCTTGGCATTCCGGTACTGTTGGTAAAGCCACC
    Luciferase
ATGGAAGACGCCAAAAACATAAAGAAAGGCCCGCGCCATTCTATCCGCTGGAAGATGGAACCGCTGGAGA
GCAACTGCATA

```

Shelf Life

The shelf life of this product is 12 months.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com