



Signosis


Innovative Plate Assay Solutions

Product Information

pE47-Luc Reporter Vector (LR-2032)

Product

pE47-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of E47 in cells.



```

                                     NheI
TTTTCTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGC-
      E47 consensus sequence
CCGGCAGGTGTCCCCCGGCAGGTGTCCCCCGGCAGGTGTCCCAGATCTGGGTATATAATGGAAGCTTG
                                     BglII                               HindIII
                                     Luciferase
GCATTCCGGTACTGTTGGTAAAGCCACCATGGAAGACGCCAAAAACATAAAGAAAGGCCCGCGCCATTCT
ATCCGCTGGAAGATGGAACCGCTGGAGAGCAACTGCATA
  
```

Shelf Life

The shelf life of this product is 12 months.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com