



Signosis


Innovative Plate Assay Solutions

Product Information

pEGR-Luc Reporter Vector (LR-2033)

Product

pEGR-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of EGR in cells.



```

                                NheI
TTTTCTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGC
EGR cis element                                BglII
GGATCCAGCGGGGGCGAGCGGGGGCCAGGATCCAGCGGGGGCGAGCGGGGGCCAAGATCTGGGTATATAAT
HindIII
GGAAGCTTGGCATTCCGGTACTGTTGGTAAAGC
                                Luciferase
CACCATGGAAGACGCCAAAAACATAAAGAAAGGCCCGGCCATTCTATCCGCTGGAAGATGGAACCGC
TGGAGAGCAACTGCATA
  
```

Shelf Life

The shelf life of this product is 6 months. Use it before expiration.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com