



Signosis


Innovative Plate Assay Solutions

Product Information

pFKHR-Luc Reporter Vector (LR-2037)

Product

pFKHR-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of FKHR in cells.



```

                                     NheI           FKHR Cis-element
TTTTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGC CAAAACAACAAAACAACAAAACAACAAAACA
      BglII           HindIII
AAGATCTGGGTATATAATGGAAGCTTGGCATTCCGGTACTGTT
                                     Luciferase
GGTAAAGCCACC ATGGAAGACGCCAAAAACATAAAGAAAGGCCCGGCGCCATTCTATCCGCT
GGAAGATGGAACCGCTGGAGAGCAACTGCATA
  
```

Shelf Life

The shelf life of this product is 12 months. Use it before expiration.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com