



# Signosis


Innovative Plate Assay Solutions

## Product Information

### **pGATA-Luc Reporter Vector (LR-2103)**

#### Product

pGATA-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of GATA in cells.



NheI      GATA cis element

GAACATCTCTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGC**CACTTGATAACAGAAAGTGATAAT**

BglII

**GATAACAGAAAGTGATAACTCTAGATCTGGGTATATAATGGAAGCTTGGCATTCCGGTACTGTTGGTAAAG**

CCACCATGGAAGACGCCAAAAACATAAAGAAAGGCCCGCGCCATTCTATCCGCTGGAAG

Luciferase

ATGGAACCGCTGGAGAGCAACTGCATAA

#### Shelf Life

The shelf life of this product is 6 months. Use it before expiration.

#### Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at [support@signosisinc.com](mailto:support@signosisinc.com)