



# Signosis


Innovative Plate Assay Solutions

## Product Information

### **pFREAC-Luc Reporter Vector (LR-2038)**

#### Product

pFREAC-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of FREAC in cells.



```

                                NheI       FREAC Cis Element
TTTTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGCGGTAAGTAAATAAACAGGTAAGTAAATAAAC
                                BglII       HindIII
AGGTAAGTAAATAAACAGATCTGGGTATATAATGGAAGCTTGGCATTCCGGTACTGTTGG
                                Luciferase
TAAAGCCACC ATGGAAGACGCCAAAAACATAAAGAAAGGCCCGGCCATTCTATCCGCT
GGAAGATGGAACCGCTGGAGAGCAACTGCATA
  
```

#### Shelf Life

The shelf life of this product is 12 months. Use it before expiration.

#### Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at [support@signosisinc.com](mailto:support@signosisinc.com)