



Signosis


Innovative Plate Assay Solutions

Product Information

pHIF-Luc Reporter Vector (LR-2002)

Product

pHIF1-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of HIF1 in cells.



```

                                NheI           HIF1 cis element
TTTTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGCGTGACTACGTGCTGCCTAGGTGAC
                                BglIII
TACGTGCTGCCTAGGTGACTACGTGCTGCCTAGGTGACTACGTGCTGCCTAGAGATCTGGTATATAAT

Hind III                               Luciferase gene
GGAAGCTTGGCATTCCGGTACTGTTGGTAAAGCCACCATGGAAGACGCCAAAAACATAAAGAAAGGCCCGG
CGCCATTCTATCCGCTGGAAGATGGAACCGCTGGAGAGCAACTGCATA

```

Shelf Life

The shelf life of this product is 12 months. Use it before expiration.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com