



Signosis


Innovative Plate Assay Solutions

Product Information

pHNF1a-Luc Reporter Vector (LR-2039)

Product

pHNF1a-Luc reporter vector is a luciferase-based reporter construct for quantitative measurement of activity of HNF1a in cells.



```

                                NheI           HNF1a Cis-element
TTTTCTATCGATAGGTACCGAGCTCTTACGCGTGCTAGC GCCCGTGGTTAGTAATTACTAAGGTTAGTAA
                                BglII           HindIII
TTACTAAGGTTAGTAATTACTAAAGATCTGGGTATATAATGGAAGCTTGGCATTCCGGTACTGTT
                                Luciferase
GGTAAAGCCACC ATGGAAGACGCCAAAAACATAAAGAAAGGCCCGGCGCCATTCTATCCGCT
GGAAGATGGAACCGCTGGAGAGCAACTGCATA
  
```

Shelf Life

The shelf life of this product is 12 months. Use it before expiration.

Contacting Signosis

For technical questions, contact our technical support group by telephone at 1-408-747-0771 or by email at support@signosisinc.com